

PURTAB FOGGER

The YMCA has a hospital grade disinfectant that we will be using daily. It is mixed with water and then put into a sprayer. The product turns into a mist and will be sprayed onto all of the workout equipment. The product must remain on the surface for 1 minute and then air dry.

PURTABS are effervescent disinfecting and sanitizing tablets that deliver high performance in a safe and effective formula. PURTABS dissolve in water to produce a variety of strengths that meet the range of needs for every facility. From the low-level sanitization required for food-contact surfaces, to powerful infection control for C. diff in health care settings.

The active agent in PURTABS is sodium troclosene (NaDCC). While it does provide chlorine in the form of hypochlorous acid (HOCL) for sanitizing and disinfection purposes, it is not a hypochlorite like chlorine bleach. The chemistry and mode of action of NaDCC is significantly different producing a solution that is stable once diluted, particularly in the presence of organic contaminants.